
## Wondfo® 2019-nCoV Antigentest (Lateral-Flow-Methode)

### Nasen-Rachen-Abstrich / Rachenabstrich

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

Rev A1  
Ref. 2021.01.18

**Contents:**

1. 1 Individual sealed pouches, each pouch contains:
  - 1 x Test cassette
  - 1 x Desiccant pouch
2. 1 Extraction Buffer (400  $\mu$ L/tube)
3. 1 Sample Extraction Tube
4. 1 Dropper
5. 1 Nasopharyngeal Swab **STERILE**
6. 1 Procedure Card
7. Instructions for use

### Wondfo® Leading POCT Manufacturer

## 2019-nCoV Antigen Test (Lateral Flow Method)

LOT

REF

WI96M163

WI96P0005

2021-11-26

2020-11-27

FOR IN VITRO DIAGNOSTIC USE ONLY

# PROCEDURE

## 2019-nCoV Antigen Test

### INTENDED USE

Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) is an immunochromatographic for rapid, qualitative detection of novel coronaviruses (2019-nCoV) antigen extracted from nasopharyngeal swab or oropharyngeal swab specimen. The test is to be used as an

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CERTIFICATE OF CONFORMITY

新冠病毒抗原检测试剂盒 (胶体金法)  
Wondfo 2019-nCoV Antigen Test (Lateral Flow Method)

产品名称 Product Name	新冠病毒抗原检测试剂盒 (胶体金法) Wondfo 2019-nCoV Antigen Test (Lateral Flow Method)
产品货号 Product Code	WI96
生产批号 Batch Number	WI96M163
生产日期 Manufacture Date	2021-11-26
有效期 Expiry Date	2020-11-27
执行标准 Specification	YY 0451-2018
检测方法 Detection Method	胶体金法 Gold Particle Test
生产厂家 Manufacturer	佛山沃丰生物科技有限公司 Foshan Wondfo Biotech Co., Ltd.
厂家地址 Address	广州市番禺区科学城珠村山脚8号 No.8 Lushan Road, Science City, Luogang District, 510663, Guangzhou, P.R. China

EC DECLARATION OF CONFORMITY  
According to the In Vitro Diagnostic Medical Device Directive 98/79/EC

<b>Manufacturer:</b>	Guangzhou Wondfo Biotech Co. Ltd.	
<b>Address:</b>	No.8, Lizhishan Road, Science City, Luogang District, 510663, Guangzhou, P.R. China	
<b>In vitro diagnostic device(s):</b>	<b>Product Name:</b>	<b>Cat. No.:</b>
	Wondfo 2019-nCoV Antigen Test (Lateral Flow Method)	W196
	<b>IVDD Classification:</b>	Other, for professional use
This declaration of conformity is issued under the sole responsibility of the manufacturer that that the above product(s) meet(s) the provisions of the European Directive 98/79/EC for in vitro Diagnostic Medical Devices.		
The following (harmonized) standards have been applied:		
EN ISO 13485: 2016	EN ISO 14971: 2012	EN 13612:2002
EN ISO 15223-1:2016	EN ISO 18113-1: 2011	EN ISO 18113-2: 2011
EN ISO 23640: 2015	EN 13641: 2002	EN 62366: 2008
The conformity with the requirements of the Directive has been assessed following the procedure(s) outlined in the following annexes of the Directive: <u>Annex III, excluding 6</u>		
<b>Notified Body (if consulted):</b>	Not applicable.	
Technical documentation demonstrating compliance is kept by the manufacturer and can be made available by the authorized representative in Europe:		
Qarad BV, Ciplastraat 3, 2440 GEEL, Belgium		
<i>Guangzhou Nov. 6, 2020</i>	Yaqin Chi, Regulatory Affairs Director <i>Yaqin Chi</i>	
(Place and date of issue)	(name and signature or equivalent marking of authorized person)	

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Guangzhou Wondfo Biotech Co., Ltd.